

**A.R.C.
LASER**

enlighten your surgery.

FOX

BASIC SURGICAL LASER

ENT
PERFECTION



CLINICAL
USE
OT

IN OFFICE
USE
OPD

**980
nm**

POWERFUL

ALLROUNDER

**MINIMAL
INVASIVE**




FOX

Cuts without almost
any bleeding

Rapidly coagulates tissue

Should be present in
any physicians office

Is available at a
fraction of the cost
of other lasers



HiFlex fibers for max freedom



Light weight flexibility



Noiseless system

Fiber guided freedom and flexibility

The lightweight and portable FOX laser fits in the palm of your hand. The ultra powerful setup allows laser beams with up to 12 W – perfect to perform surgery or treatments in outpatient setup.

Order Information
FOX 980 ENT

Best-in-Class Versatility

A.R.C. Laser is the synonym for laser technology "Made in Germany". All research and development for our lasers & optical equipment takes place in Germany.

SPECIFICATIONS

	FOX 980
Laser Wavelength	980 nm
Output Power (Laser)	12 W
Operating Mode	Pulsed, Single: 0.1 ms to CW
Beam Delivery	Fiber @ UC connector A.R.C.
Display / Control	Color touch screen / Foot switch
Cooling	Internal, forced air
Aiming Beam	Green 532 nm, adjustable, <1mW
Power Requirement	100-240 V AC, 47/63 Hz, 1.06-0.45 A integrated batteries
Dimension (HxWxD)	17.4 cm x 14.2 cm x 16.3 cm
Weight (Laser)	1.2 kg

**A.R.C.
LASER**

LASER...INNOVATION
MADE IN GERMANY

A.R.C. Laser GmbH
Bessemersstr. 14
90411 Nuremberg

+49 911 217 79 -0
+49 911 217 79 99
info@arclaser.com

distributed by:

InPhase Medical